


Sling Solutions

Seamless integration


Hygiene Sling




General Purpose Sling



Positioning Sling



Back Belt Sling




Full Body Support Sling

BARTON™
MEDICAL CORPORATION


General Slings

General slings are multi-purpose strap-and-loop slings with built-in flexibility to achieve leg and body positioning. General slings have equal application on ceiling lifts and floor lifts. The non-proprietary strap-and-loop attachments are industry standard and suitable for use on any carrybar designed for strap-and-loop slings.


Straps and Loops




General Purpose




Model	Style	Size
NSB-100-LM	Mesh	L
NSB-100-MM	Mesh	M
NSB-100-SM	Mesh	S
NSB-100-LP	Padded	L
NSB-100-MP	Padded	M
NSB-100-SP	Padded	S




General Purpose with Headsupport




Model	Style	Size
NSB-100-LMH	Mesh	L
NSB-100-MMH	Mesh	M
NSB-100-SMH	Mesh	S
NSB-100-LPH	Padded	L
NSB-100-MPH	Padded	M
NSB-100-SPH	Padded	S




Full Body Support




Model	Style	Size
NSB-200-LM	Mesh	L
NSB-200-MM	Mesh	M
NSB-200-SM	Mesh	S
NSB-200-LP	Padded	L
NSB-200-MP	Padded	M
NSB-200-SP	Padded	S



Hygiene



Model	Style	Size
NSB-300-LB	Buckle	L
NSB-300-MB	Buckle	M
NSB-300-SB	Buckle	S
NSB-300-LV	Velcro	L
NSB-300-MV	Velcro	M
NSB-300-SV	Velcro	S




Sling Dimensions

	General Purpose			General Purpose with Headsupport			Full Body Support			Hygiene		
	S	M	L	S	M	L	S	M	L	S	M	L
Width	32½"	38½"	44½"	32½"	38½"	44½"	28"	34"	40"	32"	38"	44"
Commode	10"	18"	21"	10"	18"	21"	6½"	10½"	12½"	15½"	21½"	27½"
Length	42½"	45½"	50"	53½"	58½"	65"	49½"	52½"	56½"	37"	38"	39"
Body	22½"	24"	26½"	33½"	37"	41½"	35½"	37½"	39½"	10"	10"	10"

Extension Loops

Extension loops thread through permanent loops to allow secure field modifications to the length of any strap. Extension loops are also an economical alternative to a secondary carrybar on a portable ceiling lift.



Special Slings

Positioning Sling


The sling is used for re-positioning and turning and for supine transfers. In terms of public awareness, caregiver injuries from patient handling are commonly associated with lift-and-transfer activities, but institutional experience indicates that as much as 50 percent of injuries may be caused by re-positioning and turning.

Positioning slings are used primarily with ceiling lifts, but the standard sling may also be used with a floor lift for lifting a fallen resident.


The three models are universal fit, but serve different purposes:

- The dressing sling is a short positioning sling used for turning and dressing.
- The standard sling is a general purpose positioning sling.
- The extra-long sling is for multiple repositioning of restless residents.

Positioning Sling




Model	Style
NSB-400-D	Dressing
NSB-400-S	Standard
NSB-400-XL	Extra long




Back Belt

The back belt is a specialty sling for use on a sit-to-stand device.



Model	Style
NSB-500-B	Buckle
NSB-500-V	Velcro







Custom Slings

All standard slings can be customized for *individual elements* including commode opening, strap length, the number of loops, the dimensions of the body of the sling, the location and thickness of padding, the color of binding, the number and placement of grab handles, seatbelt or Velcro attachments.


Secondly, standard slings can be customized for *overall size* from XS to XXL.

Thirdly, slings can be made of *alternate materials* for comfort, infection control or other purposes.


In addition, slings can be *designed from scratch* to meet individual requirements. Examples of commonly requested custom slings are bathing slings, stretcher slings, bath chair slings, walking slings and operating room slings.



Positioning the legs - Full Body Support




Legs divided -
cross-over style




Legs relaxed divided -
center hook configuration

Positioning the body - General Purpose




Upright position -
“Legs long shoulders
short”



Reclining position -
“Legs short shoulders long”

Standard sling sizes



XS, XL and XXL up to 600 lbs are
available custom.

Sling binding colors

Small	purple
Medium	grey
Large	blue
Extra large	black

Sling loop colors

Small	purple
Medium	grey
Large	black

Binding and loops are color-coded to facilitate size selection and level loading.

Sling materials

	Body	Straps	Padding
Flame retardant polyester	X		
High-strength nylon webbing		X	
Closed-cell polyethylene			X

Barton Medical slings do not contain latex.

Sling-to-user compatibility

Sling Type	User ability				
	Able to stand	Upper body	Dependent	Highly involved	Most users
General Purpose Sling			X		
Full Body Support Sling				X	
Hygiene Sling		X			
Positioning Sling					X
Back Belt	X				

Sling-to-equipment compatibility

Sling Type	Lift type		
	Floor lift	Ceiling lift	Sit-to-stand
General Purpose	X	X	
Full Body Support	X	X	
Hygiene	X	X	
Positioning		X	
Back Belt			X

Distributed by:



The information in this brochure is intended as a general guideline. Sling selection in each individual instance should be done in consultation with competent medical authority. Refer to Barton Medical's Sling Owner's Manuals for complete instructions on the use and the care of slings.

Barton Medical Corporation 5725 Hwy 290 West Suite #103 Austin, TX 78735

Toll free (877)8-BARTON or 512-476-7199 Fax 512-476-7190

E-mail sales@BartonMedical.com

Website www.BartonMedical.com