



CONVERTIBLE™

H-Series



H-250 (lateral transfer position)



H-250 with scale option



H-250

BARTON™
MEDICAL CORPORATION
www.bartonmedical.com

“Your Safest Move.”

Tel: 1-800-779-5611

Fax: 1-866-669-3702

The Barton *Convertible™ H-250* multi-positional chair helps provide a no-lift solution for patient transfers in a safe, dignified environment for both patient and caregiver. Easily converting between a stretcher and an upright mobile chair in seconds, the Barton *Convertible™ H-250* enables caregivers to comfortably transfer and transport a patient weighing up to 250 pounds. The *Convertible™ H-250* adapts from stretcher to transport chair with an infinite number of positions available, including recline, Tilt-in-Space and Trendelenburg capabilities. Design innovations include a fully adjustable headrest, heavy-duty steel frame, and optional scale to make weighing mobility-impaired patients easy, safe, and comfortable.



H-250



CONVERTIBLE™

H-Series



H-250 with Scale Option

Dramatic advances in medical technologies have been made in recent years. However, the difficult and awkward process of transferring a patient from bed to chair and back again is still a strenuous and dangerous process that frequently results in injury to both patient and caregiver.

Barton® Convertible™ H-250

The Barton® Convertible™ H-250 chair and Positioning and Transfer System™ (PTS™) provide the “No-Lift” solution that is revolutionizing healthcare by providing a safe, comfortable and dignified method of transferring and moving non- and limited-ambulatory patients. Tested and proven in hospitals, rehab facilities, nursing homes and home-care environments, the Barton® Convertible™ H-250 can help reduce both time and caregivers required, while significantly reducing injuries and costs.


Standard Features

- Reclines to stretcher position to transfer between bed and chair
- Tilt-in-Space capability, providing pressure relief
- Fully adjustable headrest with side-support
- Unlimited reclining positions available for patient comfort
- Trendelenburg position capability
- Sturdy steel frame accommodates patients weighing up to 250 lbs.
- Overall 24-inch chair width can move through standard doorways

Options

- Integrated digital scale conveniently weighs patients in seconds, accurate to 1/10 lb.
- IV Pole
- O₂ Tank Holder

- └ upright position
- ✓ reclined position
- stretcher position
- Trendelenburg position



The Barton® Convertible™ H-Series easily adjusts to a variety of comfortable positions: upright, reclined, stretcher, and Trendelenburg position.



A Fully adjustable headrest includes side support.



B Comfortable Tilt-in-Space positioning.




C Sturdy 6-inch casters maneuver gently and easily on carpet, tile or hardwood floors.


The Barton® Positioning and Transfer System™




H-250 with accompanying Barton® PTS™.




1 Place chair at bed side in horizontal position.



2 Attach sheet to the Barton® PTS™.



3 Transfer patient by turning handle.



4 Incline chair and you're ready to go!

Process showing patient transfer. To transfer from the chair to the bed, reverse the process.

Call Barton® or your local representative for information.

Tel: 1-800-779-5611

Fax: 1-866-669-3702

BARTON
MEDICAL CORPORATION
www.bartonmedical.com

