



## Users “see the light” with the KODAK Miniloader 2000P

*Automated roomlight film handling for mammography  
Means more time with patients, reduced artifacts, fewer retakes*

*“With the Kodak miniloader 2000P, we can handle a higher volume of patients and preserve image quality. Artifacts are very rare at our facility.*

*“We image almost 40 patients a day with just two exam rooms. That’s not possible with a darkroom, but it is possible with a Kodak miniloader 2000P.*

*“The Kodak miniloader 2000P offers excellent productivity, image quality, and reliability. That’s why we ordered a second one.”*

**Debbie Mannella**, RT (R) (M)

*Mammography Technologist*

Marymount Hospital,  
Garfield Heights, Ohio

*"The Kodak miniloader 2000P is convenient and easy to use. It allows me more time to spend with patients...which is much more productive than standing in a darkroom."*

*"The hands-free operation eliminates fingerprints and other artifacts, which reduces the need for retakes."*

**Joyce Bemiller**, RT (R) (M)

*Mammography Technologist*

Samaritan Women's Health Services,  
Ashland, Ohio

*"The Kodak miniloader 2000P enables us to have a processing area close to the mammography suite where there is no darkroom. That's very convenient."*

*"The miniloader enhances both patient throughput and patient care. We spend time talking to patients instead of handling film."*

**Alice Jackson**

*Director of Medical Imaging*

Salem Community Hospital,  
Salem, Ohio

## SPECIFICATIONS

### Cycle time

12 to 15 sec., depending on film size and magazine position

### Dimensions

2000P Configuration with integrated processor:

Width: 73 cm/28.3 in.

Depth: 78 cm/31.5 in.

Height: 138 cm/55 in.

Weight: 250 kg/550 lbs.

### Electrical

110/240 V, single phase

300 W max. consumption

UL listed; CSA, TUV-GS, and CE certified

### Cassettes

Accepts Kodak Min-R 2 cassettes (18 x 24 cm and 24 x 30 cm). Can be modified for Kodak X-Omatic video film holders (8 x 10 in.).

### Accessories

- Kodak miniloader 2000P humidifier
- Kodak miniloader 2000P/Min-R processor ready kit

Note: Processor and accessories are sold separately.

Health Imaging Division  
EASTMAN KODAK COMPANY  
Rochester, NY 14650

KODAK CANADA INC.  
Toronto, Ontario M6M 1V3  
CANADA

Outside the U.S. or Canada, please contact your local Kodak company.

[www.kodak.com/go/mammo](http://www.kodak.com/go/mammo)

