



Auf diesen Seiten weisen wir auf wichtige aktuelle Beiträge in unseren Schwesterzeitschriften hin. Wenn Sie die Seiten online lesen, dann können Sie

die Artikel mit einem Klick direkt aufrufen, ansonsten sind sie durch Eingabe der DOIs über Wiley Online Library leicht online zugänglich.




Cross-Coupling Reactions

B. Elsler, A. Wiebe, D. Schollmeyer, K. M. Dyballa, R. Franke, S. R. Waldvogel*

Source of Selectivity in Oxidative Cross-Coupling of Aryls by Solvent Effect of 1,1,1,3,3,3-Hexafluoropropan-2-ol

A remarkable and tunable solvent effect is responsible for the highly selective cross-coupling of phenols with arenes, phenols or aniline derivatives (see scheme). The solvent mixture allows the decoupling of nucleophilicity and oxidation potential in some extent.



Chem. Eur. J.
DOI: [10.1002/chem.201501604](https://doi.org/10.1002/chem.201501604)




Trifluoromethylation


H. Egami, Y. Usui, S. Kawamura, S. Nagashima, M. Sodeoka*

Product Control in Alkene Trifluoromethylation: Hydrotrifluoromethylation, Vinylic Trifluoromethylation, and Iodotrifluoromethylation using Togni Reagent

Three for the price of one: Hydrotrifluoromethylation, vinylic trifluoromethylation, and iodotrifluoromethylation of simple alkenes have been selectively achieved by using Togni reagent in the absence of any transition metal catalyst with the appropriate salts and solvents.



Chem. Asian J.
DOI: [10.1002/asia.201500359](https://doi.org/10.1002/asia.201500359)




Molecular Probes


S. Hou, S. E. Johnson, M. Zhao*

A One-Step Staining Probe for Phosphatidylethanolamine

Phospholipid probe: Cellular characterization of the phospholipid phosphatidylethanolamine (PE) is challenging owing to the absence of proper molecular tools. We constructed a one-step molecular probe based on the lantibiotic duramycin that retained the PE binding specificity but not the cytotoxicity of the native compound, thus making it an attractive tool for future PE characterization studies.



ChemBioChem
DOI: [10.1002/cbic.201500127](https://doi.org/10.1002/cbic.201500127)



ChemPhysChem



DOI: 10.1002/cphc.201500321

Li-Ion Chemistry

D. Wohlmuth,* V. Epp, M. Wilkening

Fast Li Ion Dynamics in the Solid Electrolyte $\text{Li}_7\text{P}_3\text{S}_{11}$ as Probed by ^6Li NMR Spin-Lattice Relaxation

Extremely mobile: Li ions in the glass ceramic $\text{Li}_7\text{P}_3\text{S}_{11}$ are distributed over seven crystallographically inequivalent sites. ^7Li NMR relaxometry is used to reveal the resulting complex Li self-diffusivity from a microscopic point of view. Heterogeneous ion dynamics is characterized by a broad distribution of jump rates and activation energies.



ChemMedChem


DOI: 10.1002/cmdc.201500163

Prodrugs

M. Zindler, B. Pinchuk, C. Renn, R. Horbert, A. Döbber, C. Peifer*

Design, Synthesis, and Characterization of a Photoactivatable Caged Prodrug of Imatinib

Don't fence me in! A caged prodrug of imatinib, the first approved kinase inhibitor, was designed, synthesized, and characterized for photokinetic properties. By using a PDGF-Rβ assay, we demonstrated that upon UV irradiation at $\lambda = 365 \text{ nm}$ the biological activity of imatinib is restored.



ChemSusChem



DOI: 10.1002/cssc.201500638

Lithium Storage

K. Kim, J. Moon, J. Lee, J.-S. Yu, M. Cho, K. Cho, M.-S. Park,*
J.-H. Kim,* Y.-J. Kim

Mechanochemically Reduced SiO_2 by Ti Incorporation as Lithium Storage Materials

Taking charge! A mechanochemically reduced $\text{SiO}_{2-x}/\text{TiO}_{2-x}$ composite is designed by thermodynamic calculations and synthesized by high-energy mechanical milling as a promising anode material for lithium-ion batteries. By incorporation of Ti, amorphous SiO_2 can be effectively reduced during the synthesis, leading to the in situ formation of a $\text{SiO}_{2-x}/\text{TiO}_x$ composite that shows a high reversible capacity with excellent cycle performance.



ChemCatChem


DOI: 10.1002/cctc.201500192


Electrocatalysis

M. Seredych, K. László, T. J. Bandoz*

Sulfur-Doped Carbon Aerogel as a Metal-Free Oxygen Reduction Catalyst

Pores for thought: The thermal treatment of a carbon aerogel with H_2S results in the introduction of thiophenic compounds to the carbon matrix, which increased the efficiency of the oxygen reduction reaction. As a result of their hydrophobicity, the aerogels withdraw O_2 from the electrolyte. The specific micro-/mesoporosity enhances the accessibility of the surface sites to oxygen dissolved in water.






Surfactants


N. Barbero, C. Magistris, P. Quagliotto,* L. Bonandini, C. Barolo, R. Buscaino, C. Compari, L. Contardi, E. Fisicaro, G. Viscardi

Synthesis, Physicochemical Characterization, and Interaction with DNA of Long-Alkyl-Chain Gemini Pyridinium Surfactants

Bringing it all together: This study on gemini pyridinium surfactants (see figure) contributes to the knowledge about the aggregation behavior of these compounds and their interactions with biological molecules. Structure–activity relationships were found in the context of the physico-chemical properties of these compounds and DNA compaction. These insights should improve the understanding of surfactant-based gene delivery.



ChemPlusChem
DOI: 10.1002/cplu.201500007




Molecular Wires

A. Mulas, Y.-M. Hervault, L. Norel, S. Rigaut,* C. Lagrost*

Electron-Transfer Kinetics in Polymetallic Carbon-Rich Ruthenium(II) Bis(σ-arylacetylides) Wires Connected to Gold

Well connected: Self-assembly of organometallic molecular wires consisting of mono-, bi-, and trinuclear ruthenium bis(σ-arylacetylides) functionalized with different anchoring groups all display fast electron-transfer dynamics (10^4 s^{-1}) together with discrete oxidation events at remarkably low potentials.



ChemElectroChem
DOI: 10.1002/celec.201500206




Cooperative Phenomena


M. Wang, T. Weyhermüller, K. Wieghardt*

Determining the Electronic Structure of a Series of $[(\text{phen})_3\text{M}]^0$ ($\text{M} = \text{Ti}, \text{V}, \text{Mo}$) and $[(\text{pdi})_2\text{M}]^{n+}$ ($\text{M} = \text{Cr}, \text{Mo}$) Complexes: Coordination of Neutral Ligands vs. π -Radical Anions

The appropriate oxidation levels of the tridentate ligands pyridine-2,6-diimine (pdi) and 2,9-dimethyl-1,10-phenanthroline ($^{\text{Me}}\text{phen}$) as neutral species, monoanion π -radical and dianion in “low valent” octahedral complexes $[\text{M}(^{\text{Me}}\text{phen})_3]^0$ ($\text{M} = \text{Ti}, \text{V}, \text{Cr}, \text{Mo}$) and $[\text{Mo}(\text{pdi})_2]^n$ ($n = 0, 1, +$) have been experimentally determined by X-ray crystallography and by UV/Vis and EPR spectroscopy.



Eur. J. Inorg. Chem.
DOI: 10.1002/ejic.201500420




Glyconanodiamonds


C. Fessele, S. Wachtler, V. Chandrasekaran, C. Stiller, T. K. Lindhorst,* A. Krueger*

Thiourea-Bridged Nanodiamond Glycoconjugates as Inhibitors of Bacterial Adhesion

Thiourea-bridging is a new approach to conjugating nanodiamond nanoparticles and biomolecules. This was exemplified by preparation of saccharide-functionalized nanodiamonds (NDs) that were characterized and employed in biological studies. These studies revealed specific ND binding to *E. coli* bacterial cells.



Eur. J. Org. Chem.
DOI: 10.1002/ejoc.201500532




ChemistryOpen
DOI: 10.1002/open.201500053

Chemosensors

F. Sancenón,* L. Pascual, M. Oroval, E. Aznar, R. Martínez-Máñez*

Gated Silica Mesoporous Materials in Sensing Applications

Heightened senses! Capped silica mesoporous supports (SMSs) have found wide application as delivery vehicles due to their ability to release their cargo in response to chemical, biochemical or physical stimuli. Here, a review of the application of gated SMSs in recognition protocols for anions, cations, small molecules, and large biomolecules is presented.





Asian J. Org. Chem.
DOI: 10.1002/ajoc.201500248

Rotaxanes

J. Nishiyama, Y. Makita, N. Kihara*

Rapid and Efficient Acylative Active Transport on a Rotaxane

Very moving: Smooth and efficient active transport was achieved on a rotaxane consisting of a secondary ammonium salt bearing a cyclopentyl group and dibenzo-24-crown-8-ether (DB24C8) mover by applying rapid acylation conditions. Thus, a two-step, unidirectional transport system was developed. The DB24C8 mover enters the track over the left cyclopentyl group. The rapid acylation actively transports the mover to the right. Finally, the DB24C8 mover dethreads over the right terminal group.




ChemNanoMat
DOI: 10.1002/cnma.201500045

Electrode Materials

V. Kumar, X. Wang, P. S. Lee*

Oxygen-Ions-Mediated Pseudocapacitive Charge Storage in Molybdenum Trioxide Nanobelts

Surface species: The coordination of surface atoms can tune the electrochemical properties of electrode materials. The oxygen species present on the surface of MoO₃ nanobelts affect their charge-storage properties.



ChemViews magazine
DOI: 10.1002/chemv.201400081

Science Communication

R. Threlfall

Presentation Tips: Making it Memorable

Making scientific presentations engaging for the audience can be a challenge. Richard Threlfall, *Asian Journal of Organic Chemistry*, gives tips on how to design your slide show, step into the spotlight, and add some magic to your presentation.

