

workflow. Neutral loss scanning, including multistage activation, on an LTQ-orbitrap mass spectrometer is an alternative to the precursor ion scanning used on triple quadrupole mass spectrometers. This high resolution and high mass accuracy mass spectrometer can be used successfully for the identification of protein phosphorylation sites and the identification sites in the protein UCP1 will be used as an example.

doi:10.1016/j.bbabo.2008.05.417

S16/3 Mammalian liver cytochrome *c* is tyrosine-48 phosphorylated *in vivo*, inhibiting mitochondrial respiration

Hong Yu^a, Icksoo Lee^a, Arthur R. Salomon^b, Kebing Yu^b, Maik Hüttemann^a

^aWayne State University School of Medicine, Detroit, Michigan, USA

^bBrown University, Providence, Rhode Island, USA

E-mail: mhuttema@med.wayne.edu

Cytochrome *c* (Cyt *c*) is part of the mitochondrial electron transport chain (ETC), accepting electrons from *bc*₁ complex and transferring them to cytochrome *c* oxidase (CcO). Considering its central role in life (respiration) and death (apoptosis) decisions one would expect tight regulation of Cyt *c* function. We have recently shown that Cyt *c* isolated from cow heart tissue is phosphorylated on tyrosine 97 *in vivo*, which leads to inhibition of respiration in the reaction with CcO. In this study we isolated Cyt *c* from a different organ, cow liver, under conditions preserving the physiological phosphorylation state. Western analysis with a phospho-tyrosine specific antibody suggested that liver Cyt *c* is phosphorylated. Surprisingly, the phosphorylation site was unambiguously assigned to Tyr-48 by immobilized metal affinity chromatography/nano-liquid chromatography/electrospray ionization mass spectrometry, and not to the previously identified phospho-Tyr-97 in cow heart. As is true of Tyr-97, Tyr-48 is conserved in eukaryotes. As one possible consequence of Tyr-48 phosphorylation we analyzed the *in vitro* reaction kinetics with isolated cow liver CcO revealing striking differences. Maximal turnover of Tyr-48 phosphorylated Cyt *c* was 3.7 s⁻¹ whereas dephosphorylation resulted in a 2.2 fold increase in activity to 8.2 s⁻¹. Effects of Tyr-48 phosphorylation based on the Cyt *c* crystal structure are discussed.

doi:10.1016/j.bbabo.2008.05.418

S16/4 Mitochondrial comparative proteomics: Strength and pitfalls

Francis E. Sluse

Laboratory of Bioenergetics and Cellular Physiology, University of Liege, Belgium

E-mail: F.Sluse@ulg.ac.be

In this review, we describe the various techniques available to carry out valid comparative proteomics, their advantages and their disadvantages according to the goal of research. The two-dimensional gel electrophoresis and the 2D-DIGE (Differential in-gel electrophoresis) are compared with the shotgun proteomics and SILE (Stable isotopic labeling experiments). We deliver our opinion on their best fields of application in the domain of the comparative proteomics. We underline the utility of these new tools, providing mass data to study physiology and mitochondrial plasticity in front of a specific mitochondrial insufficiency or exogenous stresses. We illustrate our matter with results obtained in our laboratory specifying the importance of an approach of comparative proteomics combined on mitochondria and the cell which makes it possible to obtain important information on the statute of the mitochondrial function at the cellular level. Finally, we draw attention to the dangers of the extrapolation

of the data of proteomics to metabolic flows which require the greatest care.

doi:10.1016/j.bbabo.2008.05.419

(S16) Proteomics and mitochondria symposium abstracts (poster and raised abstracts)

S16.5 High sensitivity identification of membrane proteins by MALDI TOF-mass spectrometry using polystyrene beads

M. Zoonens, N. Bensalem, S. Masscheleyn, J. Mozo, B. Vallée, F. Brouillard, S. Trudel, D. Ricquier, A. Edelman, I.C. Guerrero, B. Miroux
University Paris Descartes and CNRS Unit 9078, Faculty of Medicine Necker-Enfants Malades, Paris, France

E-mail: mzozenens@necker.fr

Uncoupling proteins (UCPs) are transporters present in the inner membrane of mitochondria. They are found in all mammals and in plants, and they belong to the family of anionic mitochondrial carriers. Their roles are quite well established for UCP1 in nonshivering thermogenesis process, but are still in debate for the others UCPs. We study these proteins by biophysical approaches such as mass spectroscopy (MS). Peptide mass fingerprinting methods using techniques such as MALDI-TOF MS, for example, have become an important analytical tool in the identification of proteins. However, PMF of membrane proteins is a real challenge for at least three reasons. First, membrane proteins are naturally present at low levels; second, most of the detergents strongly inhibit proteases and have deleterious effects on MALDI spectra; and third, despite the presence of detergent, membrane proteins are unstable and often aggregate. We improved and showed that differential acetonitrile extraction of tryptic peptides combined with the use of polystyrene Bio-Beads triggered high resolution of the MALDI-TOF identification of mitochondrial membrane proteins solubilized either with Triton-X100 or CHAPS detergents. The sequence coverage of UCP1 obtained by this approach is close to 90% and allows further investigations for characterizing regulators bound to UCP1.

doi:10.1016/j.bbabo.2008.05.420

S16.6 Analysis of synaptic and non-synaptic mitochondria using colorless- and blue-native PAGE

Laura Brady, Rashmi U. Pathak, Martha Motherway Gildea, Keith F. Tipton, Gavin P. Davey

School of Biochemistry and Immunology and Trinity College Institute of Neuroscience, Trinity College Dublin, Ireland

E-mail: labrady@tcd.ie

Mitochondrial dysfunction contributes to a wide range of human diseases, including neurodegenerative diseases. Proteomic technology can give novel insights into the structure and composition of the brain mitochondrial electron transport chain (ETC) complexes. Using beef heart mitochondria as a control and reference point, this study reports a comparison of rat synaptic and non-synaptic mitochondrial protein profiles using blue-native and colorless-native (CN) gel electrophoresis combined with Tricine-SDS PAGE and MALDI-TOF mass spectrometry. BN-PAGE was found to be a straightforward tool for proteomic analysis of ETC complexes (I, III, IV and V) and especially for the identification of very hydrophobic membrane protein constituents that are not accessible by common isoelectric focusing/sodium dodecyl sulphate gel electrophoresis. The introduction of CN-PAGE into a three dimensional electrophoresis greatly improved the isolation and