

The proton-pumping NADH:ubiquinone oxidoreductase (complex I), the first complex of the respiratory chains, couples the transfer of electrons from NADH to ubiquinone with the translocation of protons across the membrane. The soluble NADH dehydrogenase fragment represents the electron input part of complex I and comprises 3 subunits, one flavin mononucleotide and up to six iron–sulfur clusters. Previous electron microscopic studies revealed NADH induced large conformational changes for complex I. Structural information on the soluble fragment about these changes are not available. ATR-FTIR spectroscopic investigations of the conformational change of the NADH dehydrogenase fragment in the presence and absence of NADH substrate have been applied. The secondary structure determination was based on the analysis of the intense amide I band between 1710 to 1600 cm⁻¹ which predominantly includes the C=O stretching vibration of the protein backbone. The highly organized α helical and β sheet components are negligibly altered upon NADH addition whereas variations for turns and random elements are obtained. During perfusion-induced investigations an increase of the amide I and amide II are observed which indicates a conformational change upon NADH binding. Furthermore, investigations on the whole complex I are presented.

doi:10.1016/j.bbabo.2008.05.152

S4.22 RP-HPLC AND ESI/MS/MS subunit analysis of bovine heart complex I: Application for identifying ROS-induced modifications

Andrej Musatov, Patrizia Lemma-Gray, Susan T. Weintraub, Neal C. Robinson

Department of Biochemistry, The University of Texas Health Science Center, 7703 Floyd Curl Drive, San Antonio, TX 78229-3900, USA

E-mail: robinson@uthscsa.edu

An effective method was developed for isolation and analysis of bovine heart Complex I subunits. The method utilizes C18 reversed-phase HPLC and a water:acetonitrile gradient containing 0.1% trifluoroacetic acid. Employing this system, 36 of the 45 Complex I subunits eluted in 28 distinct chromatographic peaks. Nine subunits did not elute including B14.7, MLRQ and the seven mitochondrial-encoded subunits. The method, with UV detection, is suitable for either analytical (<50 μ g protein) or preparative (>250 μ g protein) applications. Subunit(s) eluting in each chromatographic peak were initially determined by MALDI-TOF/MS with subsequent positive identification by reversed-phase HPLC-ESI/MS/MS analysis of tryptic digests. In the latter case, subunits were identified with a 99% probability using Scaffold analysis software combined with a X! Tandem protein search engine. The RP-HPLC subunit analysis method represents a major improvement over previous separation methods with respect to resolution, simplicity and ease of application. As expected, the HPLC method, combined with mass spectrometry analysis is capable of detecting oxidatively damaged subunits after *in vitro* exposure of Complex I to peroxynitrite. The results of these experiments will help elucidate the chemical nature of oxidative damage to this key electron transport chain component. Supported by NIH GMS (24795) and Welch Foundation (AQ1481).

doi:10.1016/j.bbabo.2008.05.153

S4.23 Heterologous expression, isolation and characterization of A subcomplex of the respiratory complex I from *Aquifex aeolicus*

Markus Kohlstädt, Katerina Dörner, Ramona Labatzke, Cengiz Koç, Thorsten Friedrich

Institut für Organische Chemie und Biochemie, Albert-Ludwigs-Universität Freiburg, Germany

E-mail: markus.kohlstaedt@pluto.uni-freiburg.de

The proton-pumping NADH:ubiquinone oxidoreductase (complex I) is the first complex of the respiratory chains in many bacteria and most eukaryotes. It is the least understood respiratory complex, due to its size and intricate assembly. Complex I of the hyperthermophilic bacterium *Aquifex aeolicus* is made up of 13 different subunits named NuoA2–N1. The flavoprotein subcomplex contains subunits NuoE and F and is considered as the simplest catalytically active NADH dehydrogenase. It harbours one FMN, the binuclear iron–sulfur cluster N1a and the tetranuclear cluster N3 as well as the NADH binding site. This subcomplex was heterologously overproduced in *Escherichia coli* and isolated by standard chromatographic techniques. The proteins of the preparation were identified as NuoE and NuoF by N-terminal sequencing and mass spectrometry. The expected cofactors of the preparation were detected by EPR- and UV/vis spectroscopy. The preparation is thermostable and catalyzes NADH/ferricyanide oxidoreductase activity at temperatures up to 85 °C. Furthermore, the preparation is homogenous as shown by dynamic light scattering and meets all prerequisites for crystallization attempts.

doi:10.1016/j.bbabo.2008.05.154

S4.24 Modelling, cloning and expression of the *P. falciparum* NADH dehydrogenase

Nicholas Fisher, Stephen A. Ward, Giancarlo A. Biagini

Liverpool School of Tropical Medicine, Liverpool L3 5QA, UK

E-mail: n.e.fisher@liverpool.ac.uk

The respiratory chain of the human malaria parasite *Plasmodium falciparum* lacks a canonical protonmotive NADH:ubiquinone oxidoreductase (Complex I), containing instead a single-subunit, non-protonmotive Ndh2, similar to that found in plant mitochondria, fungi and some bacteria. As such, the *P. falciparum* Ndh2 (PfNdh2) presents itself as an attractive anti-malarial chemotherapeutic target, and we have developed a structural model and heterologous expression system for this enzyme to facilitate its physicochemical and enzymological characterisation. Structural modelling of PfNdh2 suggests the presence of two Rossman folds forming the flavin and NADH binding sites, with membrane attachment via a C-terminal amphipathic helix. Heterologous expression of PfNdh2 in *E. coli* NADH dehydrogenase knockout strain ANNO222 yields an active enzyme with Km values for NADH and decylubiquinone of 14- and 12 μ M respectively in crude membrane preparations. In addition, the recombinant PfNdh2 is strongly inhibited by 1-hydroxy-2-dodecyl-4(1H)quinolone (IC50=70 nM in crude membrane preparations) with decylubiquinone as substrate. Purification of His-tagged PfNdh2 for high throughput drug screening, stopped-flow spectrophotometry and NMR/FTIR spectroscopic investigation is currently underway.

doi:10.1016/j.bbabo.2008.05.155