Here we demonstrate a method that allows light-induced activation of voltage-gated ion channels and the concurrent imaging of membrane potential changes using voltage-sensitive dyes. This light-induced voltage clamp (LIVC) method uses photostimulation through channelrhodopsin-2 (ChR2) to activate voltage-gated ion channels. ChR2 allows light to be immediately transduced into a depolarizing ionic current, which in turn causes voltage-gated ion channels to open. In our system we coexpressed ChR2 either with the voltage-gated potassium channels hERG or hKv1.5 in cell lines and in Xenopus oocytes. In electrophysiological experiments we show that light-induced depolarization through ChR2 sufficed to activate hERG as well as hKv1.5 channels. We were further able to optically monitor the light-induced membrane de- and hyperpolarizations on a millisecond timescale with the voltage-sensitive RH421 and Annine6. The fluorescence readout reflected the dose-response relationships of the hERG blocker Terfenadine and the hKv1.5 inhibitor DPO-1 obtained from patch-clamp measurements.

LIVC represents a solely optical technology with remote activation of the target voltage-gated ion channels by the delivery of a flash of blue light and simultaneous detection of their activity employing voltage-sensitive dyes. It combines the high-throughput of optical methods with the high-content of patch clamp concerning high temporal resolution, membrane potential control and repetitive stimulation.

3694-Pos

Functional Studies of Volvox Channelrhodopsin Chimeras Matthias Prigge, Franziska Schneider, Satoshi S.P. Tsunoda,

Peter Hegemann.

Humboldt University Berlin, Berlin, Germany.

Recently, two new members of the channelrhodopsin family have been discovered in the multicellular green alga Volvox carteri and named VChR1 and VChR2. VChR2 shows a similar characteristic to the known Chlamydomonas ChR2 [1], which is successfully employed in optogenetics. It shows a maximum absorption at 470nm and can be well expressed in eukaryotic cells like HEK293T cells. In contrast, VChR1 shows the most red-shifted absorption peaking at 540nm. It has been successfully used to depolarize hippocampus cells even at 590nm [2].

Longer wavelengths would enable deeper brain stimulations and would therefore be less invasive and harmful for tissues. The use of VChR1 was hindred by low expression levels and membrane localization. In this study, we systematically exchanged helices, creating chimeras of VChR1 and VChR2. Functional chimeras show differences in expression levels, membrane localization and absorption. Several chimeras show higher expression levels and plasma-membrane localization in HEK293 cells than VChR1 while preserving absorption at 540nm.

Thereby we were able to identify key determinants causing the colour-shift from VChR1 to VChR2, located in the last three helices of these proteins. Mutagenesis of the relevant amino acids that are thought to be promising candidates for fine-tuning to longer wavelength absorption.

[1] Channelrhodopsins of Volvox carteri are photochromic proteins that are specifically expressed in somatic cells under control of light, temperature, and the sex inducer. Kianianmomeni A, Stehfest K, Nematollahi G, Hegemann P, Hallmann A. Plant Physiol. 2009 Sep;151(1):347-66.

[2] Red-shifted optogenetic excitation: a tool for fast neural control derived from Volvox carteri. Zhang F, Prigge M, Beyrière F, Tsunoda SP, Mattis J, Yizhar O, Hegemann P, Deisseroth K. Nat Neurosci. 2008 Jun;11(6):631-3.

3695-Pos

K2P1 Assembles with K2P3 or K2P9 to Form Sumo-Regulated Task Background Channels

Leigh D. Plant, Leandro Zúñiga, Sonia Olikara, Jeremy D. Marks, Steve A.N. Goldstein.

University of Chicago, Chicago, IL, USA.

TASK subunits K2P3 and K2P9 form homo- and hetero-dimeric channels in neurons with varied sensitivities to anesthetics, flavors and pH (Berg et al. J Neurosci. 2004. 24:6693-702; Bautista et al. Nat. Neurosci. 2008. 11:772-9). Here, K2P1 subunits (Rajan et al. 2005. Cell 121, 37-47) are shown to confer SUMO-regulation on TASK channels in rat cerebellar granule neurons (CGN). First, CFP (C) or YFP (Y)-tagged subunits studied in CHO cells by donor-decay Forster resonance energy transfer (FRET) confirm biochemical evidence for assembly of K2P2 and its native isoform (Thomas et al., Neuron 2008 58:859-70). Next, FRET shows association of C-K2P1 and Y-K2P1, Y-K2P3 or Y-K2P9 but not Y-K2P2 or Kv2.1. As expected from Rajan et al, FRET registers association of Y-SUMO and C-K2P1 but not C-K2P3 or C-K2P9. In contrast, Y-SUMO and C-K2P3 or C-K2P9 FRET when untagged K2P1 is co-expressed. Consistent with electrophysiological studies (EP) showing one SUMO per channel is sufficient to silence K2P channels (Plant et al. this meet-

ing), channels with two linked K2P1 subunits (WT-WT), WT and SUMO-insensitive K2P1 (WT-K274Q), WT and K2P3 (WT-K2P3) or WT and K2P9 (WT-K2P9), pass currents when membrane patches are exposed to SUMO-protease (SENP1) and silenced by SUMO1. In contrast, K2P3-K2P3, K2P9-K2P9, K2P3-K274Q or K2P9-K274Q are constitutively active and insensitive to SENP1 and SUMO1. Finally, in CGN, immunochemistry shows K2P1, K2P3, K2P9, SUMO, SUMO E1 conjugase and SUMO E2 ligase in plasma membrane; EP reveals IK_{SO} regulation by SENP1 and SUMO, and transfection with mutant subunits demonstrate assembly of K2P1 with K2P3 or K2P9 by FRET and EP.

Excitation-Contraction Coupling II

3696-Pos

The Effect of Apelin on Single Isolated Cardiac Myocytes from Wild-Type and Apelin / APJ KO Mice

Christian Bollensdorff¹, David Charo², Michael Ho², Giovanni Fajardo³, Daniel Bernstein³, Thomas Quertermous², Euan Ashley², Peter Kohl¹.

¹DPAG, Oxford University, United Kingdom, ²Division of Cardiovascular Medicine, Stanford University, CA, USA, ³Department of Pediatrics, Stanford University, CA, USA.

Apelin, an endogenously produced peptide discovered in 1998, is regarded as one of the most potent substances to increase net contractile output in isolated heart and whole animals, apparently in a load-dependent manner. Unlike other inotropic agents, apelin has not been shown to cause cardiac hypertrophy. However, little is known about its cellular effects, with only two papers on isolated cardiomyocytes, neither of which imposed external mechanical constraints. This study assesses dynamic parameters of cardiomyocytes from wild type (WT), apelin-receptor knockout (AJP-KO) and apelin knockout (Ap-KO) mice. In addition to steady-state parameters such as fractional shortening, we measured time to peak contraction (TTP), and maximum velocities of contraction and relaxation.

Apelin administration to unloaded cells did not change fractional shortening, but significantly (p<0.05) decreased TTP both in WT (93.1% \pm 2%, n=90) and Ap-KO cardiomyocytes (83.1% \pm 5%, n=14). The specificity of apelin effects was verified in AJP-KO myocytes, where no significant change in TTP was seen (98.2% \pm 2%, n=69).

After application of mechanical pre-loads, using the two-carbonfibre technique, TTP reduction was enhanced in Ap-KO myocytes ($64.6\% \pm 5\%$, n=11). To characterise apelin effects on developed tension, the ratio between the slopes of end-systolic and end-diastolic tension-length relations was obtained. This is a non-dimensional parameter, termed Frank-Starling-Gain (FSG), that is independent of cell cross-section and allows inter-individual comparisons. Under control conditions, FSG was 1.36 ± 0.05 (n=10). After application of 10nM apelin this rose to 1.96 ± 0.23 (n=6, p<0.05), indicating improved net contractile output. Interestingly, this was caused largely by a slope reduction of the end-diastolic tension-length relation. This could explain the combination of increased net force output with lack of hypertrophy induction, and suggests that apelin more appropriately be described as a positive lusitropic agent, rather than an inotropic one.

3697-Pos

Mechanism of the Spontaneous Beating of Skeletal-Based Precursor of Cardiomyocytes (SPOC)

Takayuki Miyanishi¹, Takao Shioya², Nagomi Kurebayashi³,

Takashi Murayama³, Masahiro Oono¹.

¹Nagasaki University, Nagasaki, Japan, ²Saga University, Saga, Japan,

³Juntendo University, Tokyo, Japan.

The skeletal muscle-based precursor of cardiomyocytes (SPOC), a stem cell that potentially differentiates into heart cells [Winitsky et al. (2005) PloS Biol 3(4), e87], were isolated from leg skeletal muscle of the rat. Initially the SPOC cells had round and small appearance. When cultured on the dishes coated with extracellular matrix, they differentiated into round, tear drop, or tubular shaped cells. During the course of their differentiation, the SPOC cells exhibited a rhythmic beating associated with spontaneous oscillation of the intracellular Ca²⁺. The spontaneous beating could be suppressed by Cd²⁺ and nifedipine, indicating possible contribution of the L-type Ca²⁺-channel to the development of the automaticity. The rate of the spontaneous beating was accelerated by the administration of isoproterenol. Immunocytochemistry of the SPOC cells indicated the existence of Nkx 2.5, cardiac TnT, dihydropyridine receptor, and ryanodine receptor. Action potentials during the spontaneous beating and the underlying membrane currents were analyzed by using the patch-clamp technique. Our results indicate that the SPOC cells develop their automaticity by a unique mechanism distinct from that of skeletal or heart muscle cells.