

at a 1:1 stoichiometry to test the lung transduction functionality of the peptide using the 53 kDa dendrimer as a model macromolecular cargo. Phage peptide-dendrimer conjugates were labelled with a fluorophore and characterised by ^1H NMR and quantitative amino acid analysis prior to administration into the airways of the IPRM model. The extent of absorption of PAMAM G5.5 alone equalled $17 \pm 6\%$ of lung deposited dose absorbed by 90 min. G5.5 dendrimers displaying one syn-LTP-1 peptide per polymer (termed G5.5-syn-LTP-1) displayed a 1.8-fold greater extent of absorption ($p < 0.05$) cf. G5.5 alone; G5.5 dendrimer displaying one equivalent of the syn-NB3 peptide showed no evidence of enhanced absorption ($p > 0.05$). The enhanced absorption of G5.5-syn-LTP-1 absorption was competitively inhibited by co-administration of 100-fold molar excess of syn-LTP-1 peptide ($p < 0.05$) but not by syn-NB-3 peptide ($p > 0.05$), an observation consistent with the participation of a specific receptor-mediated transport mechanism. As such the LTP-1 peptide motif may serve as a platform for enhancing macromolecule absorption from the airways.

Reference

1. Sergeeva, et al. *Adv Drug Deliv Rev* 2006;**58**:1622–54.
2. Campbell, et al. *BBRC* 1999;**262**:744–51.
3. Sakami, et al. *Pharm Res* 2006;**23**:270–9.

doi:10.1016/j.drudis.2010.09.439

A94

Differential transport of anionic PAMAM dendrimers across *in vitro* biological barriers

Ghaith Aljayyousi^{1,*}, Will Ford¹, Neil McKeown², Mark Gumbleton¹

¹ Welsh School of Pharmacy, Cardiff University, Cardiff, UK

² School of Chemistry, Cardiff University, Cardiff, UK

*Corresponding author.

E-mail: aljayyousig@cardiff.ac.uk (G. Aljayyousi).

Polyamidoamine (PAMAM) dendrimers are a class of branched polymers that have the potential to serve as drug carriers. This is primarily due to their extremely low polydispersity index, the ability to precisely control their size and charge, and the multiple functional groups that they bear on their surfaces giving the ability to conjugate a wide range of therapeutic molecules. The transport across *in vitro* biological barriers of cationic PAMAMs has been widely

studied with reports often indicating high barrier permeability, although interpretation of such data in the context of cation-induced barrier toxicity is often omitted. We are investigating the intrinsic biological activity of intact stable anionic dendrimer-drug conjugates where the dendrimer moiety not only confers a backbone for attachment of multiple pharmacological ligands but also offers a means to physically modulate *in vivo* tissue disposition, for example, affording access to intestinal submucosa but excluding BBB penetration. In this abstract we report the differential *in vitro* barrier permeability of a molecular weight series of anionic PAMAM dendrimers, that is, G1.5, 2935 Da; G3.5, 12,931 Da; G5.5, 52,907 Da which has supported our ongoing *in vivo* investigations. Dendrimers were fluorescently labelled and added to the apical surface of epithelial cell monolayers grown on a semi-permeable inserts (Transwell). Permeability coefficients (ρ) were determined for transport in the apical to basal direction. The epithelial models included the highly restrictive MDCK-I (TEER $5000 \Omega \text{ cm}^2$), the moderately restrictive Caco-2 (TEER $600 \Omega \text{ cm}^2$) and the low restrictive MDCKII (TEER $200 \Omega \text{ cm}^2$). For CACO-2 and MDCKII an inverse relationship was evident between dendrimer transepithelial transport and dendrimer molecular size, with dendrimer ρ decreasing approx. 5-fold $\text{G1.5} \Rightarrow \text{G3.5}$, and decreasing approx. 10-fold $\text{G1.5} \Rightarrow \text{G5.5}$. The permeability of the cell models to dendrimer transport declined as the paracellular restrictiveness of the monolayers increased. Indeed, for MDCKI monolayers dendrimer concentrations in the basal chamber remained at all times below the limit level of detection, but could be readily enhanced by briefly adding EDTA to the media. Nevertheless, predicted (based upon LLQ) ρ for dendrimer transport across MDCKI were at least $\times 10$ – 15 -fold lower than in the other cell models. Significantly, even for the smallest dendrimer, that is, G1.5, the maximum predicted (based on LLQ) ρ across MDCKI was no greater than 15% of the ρ obtained for the paracellular marker F-Na. Whereas ρ for G1.5 was 51% and 56% of that for F-Na in CACO-2 and MDCKII models, respectively. Biocompatibility studies show no affect of the anionic dendrimers upon overall barrier properties. The paracellular route is the major pathway of dendrimer transport across biological barriers. Stable pharmacologically active conjugates of dendrimer – drug are an interesting experimental therapeutic with potential to provide

differential tissue distribution/exclusion based upon physical characteristics.

doi:10.1016/j.drudis.2010.09.440

A95

Non-toxic, highly efficient delivery of nucleic acids into challenging cells using safectin transfection reagent

Steve Howell^{1,*}, Mark Bradley^{1,2}, Asier Unciti-Broceta^{1,2}

¹ Deliverics Ltd., 2 Blythwood Square, Glasgow G2 4AD, UK

² School of Chemistry, University of Edinburgh, West Mains Road, Edinburgh EH9 3JJ, UK

*Corresponding author.

E-mail: steve.howell@deliverics.com (S. Howell).

Deliverics Ltd. has developed a novel cationic lipid-mediated transfection reagent for DNA and siRNA delivery into both easy and challenging to transfect eukaryotic cells: SAFectin Transfection Reagent. This reagent is a water-based formulation of cationic and neutral lipids with programmed biodegradability. SAFectin allows for the highest transfection efficiency of nucleic acids into many cell types (e.g. immortalized cells, mESC, hMSC) with the simplest-to-use and fastest procedure in the market: (i) mix SAFectin and the nucleic acid (ii) followed by direct addition to cells, either in the presence or absence of serum and antibiotics. The formulation has been developed to have very low toxicity to cells and as such it is not necessary to remove or change culture medium following transfection. Combination of the highest/safest transfection rates on the market with the simplest to use protocol ensures optimal performance and fast results. The SAFectin Transfection Reagent is a universal system that outperforms competitor's products in each of the three defining features any researcher seeks in this kind of product: efficacy, safety and ease of use.

doi:10.1016/j.drudis.2010.09.441