

Linear Multinuclear Ru^{II} Photosensitizers

Simona Vaduvescu^[a] and Pierre G. Potvin^{*[a]}

Keywords: Ruthenium / Tridentate ligands / Bridging ligands / Electron transfer / Photochemistry

The preparation of 4'-p-Tolyl-2,2':6',2''-terpyridine (ttpy) on page 1767 of the original article^[1] calls for 2.9 mL of 15% aq. KOH and 29 mL of conc. NH₄OH. These quantities are reversed and should be 29 mL for KOH and 2.9 mL for NH₄OH.

[1] S. Vaduvescu, P. G. Potvin, *Eur. J. Inorg. Chem.* **2004**, 1763–1769.

Received: February 19, 2009
Published Online: March 6, 2009

[a] Department of Chemistry, York University,
4700 Keele Street, Toronto, ON M3J 1P3, Canada
Fax: +1-416-736-5936
E-mail: pspotvin@yorku.ca