

SPOTLIGHTS ...


Click Chemistry

A. Dondoni*

Triazole: the Keystone in Glycosylated Molecular Architectures Constructed by a Click Reaction

Chem. Asian J.

DOI: 10.1002/asia.200700015



Things don't fall apart: Like the wedge-shaped keystone in the central part of an arch, triazole units formed by Cu^I-catalyzed azide–alkyne coupling acts as the ligation element of glycosylated molecular fragments. This structural motif is present in various complex glycoconjugates such as glyco-clusters, glycodendrimers, glycopolymers, glycomacrocycles, and DNA strands.

Natural Products


Y. A. Elnakady, M. Rohde, F. Sasse, C. Backes, A. Keller, H.-P. Lenhof, K. J. Weissman, R. Müller*

Evidence for the Mode of Action of the Highly Cytotoxic *Streptomyces* Polyketide Kendomycin

ChemBioChem

DOI: 10.1002/cbic.200700050

The first macrocyclic polyketide proteasome inhibitor. Kendomycin is a macrocyclic polyketide recently discovered in several *Streptomyces* species. We provide evidence here that kendomycin's potent antiproliferative activity derives from inhibition of the proteasome. Kendomycin therefore represents the founding member of a new class of proteasome inhibitors with promise both in therapy and in chemical genetics.




Single-Molecule Study

B. I. de Bakker, F. de Lange, A. Cambi, J. P. Korterik, E. M. H. P. van Dijk, N. F. van Hulst, C. G. Figdor, M. F. Garcia-Parajo*

Nanoscale Organization of the Pathogen Receptor DC-SIGN Mapped by Single-Molecule High-Resolution Fluorescence Microscopy

ChemPhysChem

DOI: 10.1002/cphc.200700169



Near-field scanning optical microscopy (NSOM) was used to visualize the nanoscale spatial arrangement of the lectin DC-SIGN on the membrane of intact dendritic cells. Whereas the con-

focal image has a diffraction-limited resolution of about 350 nm (see picture, left), NSOM can visualize DC-SIGN nanodomains with a spatial resolution of less than 100 nm (right).

Glycomimetics


K. Marotte, C. Sabin, C. Préville, M. Moumé-Pymbock, M. Wimmerová, E. P. Mitchell, A. Imberty,* R. Roy*


X-ray Structures and Thermodynamics of the Interaction of PA-IIL from *Pseudomonas aeruginosa* with Disaccharide Derivatives

ChemMedChem

DOI: 10.1002/cmdc.200700100

The disaccharide α Fuc1-4GlcNAc has been used as a scaffold toward the synthesis of a series of derivatives targeting the fucose binding lectin from *Pseudomonas aeruginosa*. High affinity has been measured by titration microcalorimetry and structural studies rationalized the difference observed for the thermodynamics of binding between the natural and the synthetic ligands.





The first example of a triple-chlorido-bridged heterodimetallic dinuclear complex comprising an $(\eta^6\text{-arene})\text{Ru}^{\text{II}}$ fragment and a Group 1 metal crown ether fragment is presented. The solid-state structure was determined by X-ray crystallography. NMR investigations indicate the presence of a dynamic equilibrium in solution.


Heterodimetallic Dinuclear Complexes

C. A. Vock, R. Scopelliti, P. J. Dyson*

Synthesis and Crystallographic Characterisation of the Heterodimetallic Complex [(Dibenzo-18-crown-6)K-($\mu\text{-Cl}$)₃Ru($\eta^6\text{-p-cymene}$)]

Eur. J. Inorg. Chem.
DOI: 10.1002/ejic.200700430

Platinum favors alpha: Gase-phase reactions of size-selected platinum cluster cations with hydrocarbons reveal a preference for α -hydrogen atoms upon dihydrogen elimination (see scheme). These observations are in strong contrast to rhodium, which favors 1,2-elimination and explains why platinum is reactive towards methane whilst rhodium is not.



Heterogeneous Catalysis

C. Adlhart, E. Uggerud*

Mechanisms for the Dehydrogenation of Alkanes on Platinum: Insights Gained from the Reactivity of Gaseous Cluster Cations, Pt_n^+ $n = 1\text{--}21$

Chem. Eur. J.
DOI: 10.1002/chem.200700501



It's elemental: A water-soluble polymer bearing multiple metal-chelating ligands has been used as a tag for bioassays with inductively coupled plasma mass spectrometry. The tag was covalently conjugated to antibodies, and the polymer-antibody constructs were loaded with lanthanide ions (Ln^{3+}) and used for the simultaneous assay of five orthogonally labeled antibodies against cell surface antigens that differ in abundance by more than two orders of magnitude.

Bioassays

X. Lou, G. Zhang, I. Herrera, R. Kinach, O. Ornatsky, V. Baranov,* M. Nitz,* M. A. Winnik*

Polymer-Based Elemental Tags for Sensitive Bioassays

Angew. Chem. Int. Ed.
DOI: 10.1002/anie.200700796



On these pages, we feature a selection of the excellent work that has recently been published in our sister journals. If you are reading these pages on a computer, click on any of

the items to read the full article. Otherwise please see the DOIs for easy online access through Wiley InterScience.