







Chirasil L-Val, conditions: initial temp. 135 °C/15 min. rate 10 °C/min., final temp. 190 °C/10 min.




Chirasil L-Val, conditions: initial temp. 145 °C/15 min. rate 10 °C/min., final temp. 190 °C/10 min.



Chirasil L-Val, conditions: initial temp. 135 °C/15 min. rate 10 °C/min., final temp. 190 °C/10 min.



Chirasil L-Val, conditions: initial temp. 145 °C/15 min. rate 10 °C/min., final temp. 190 °C/10 min.



Chirasil L-Val, conditions: initial temp. 135 °C/15 min. rate 10 °C/min., final temp. 190